

## Pure-Take Pump:

- Purity protector stops fluid from reentering the bottle
- Locking Pump will not leak when engaged suitable for travel
- One-handed operation frees operator's hands for use
- Precision machined valves for a smooth consistent operation
- Cam-lock hinge keeps lid open and closed, reduces odors
- No silicone or latex perfect for ultra clean environments
- Compatible with any 38-400 threaded neck container
- Dispenses approximately 0.2cc of liquid with each pump

Pump Material	Stainless-Steel	
Locks for Travel	Yes	
Neck Size / Thread	38-400	
Lid-Type	Cam-lock hinge lid	
Gasket Material	F217 Tri-Seal Foam	

## **Pump Care:**

Some types of fluids may leave a residue in the mechanism of the dispenser. Over time, this can reduce the performance of the pump.

Occasional cleaning is recommended to maintain your unit to original efficiency.

Clean by flushing the mechanism thoroughly with warm soapy water. Flush the mechanism again with clean water to remove any soap residue. Do not boil or steam the dispenser. For long periods of inactivity, pumps should be flushed free of liquid and stored inverted in a dry place.

Request Instructions: Before using, users shall determine the suitability of the product for their intended use, and users assume all risk and liability whatsoever in connection therewith. See the Menda Terms and Conditions - Mendapump.com/Terms-And-Conditions.aspx	Specifications and procedures subject to change without notice. PURE TAKE, LOCKING, PUMP ONLY		erance is ±10%.
MENDA	3651 WALNUT AVE., CHINO, CA 91710	DRAWING	<b>DATE:</b>
	PHONE: (909) 627-2453	NUMBER	Aug.
	WEBSITE: <u>MendaPump.com</u>	35800	2016

## (43 mm)

1-11/16"

sh on lid Overfill dish

> 11/16" (17 mm)

1-1/8" (29 mm)

## Limited Warranty Warranty Evolutions, Limit of Liability and RMA